

FavorPrep™ Blood/Cultured Cells Genomic DNA Extraction Midi Kit

(For Research Use Only)

Kit Contents:

Cat. No. (preps)	FABGK002-S (2 preps)	FABGK002 (25 preps)
Proteinase K (Liquid)	450 µl	1050 µl × 4
FABG Buffer	4 ml	42 ml
W1 Buffer * (Concentrate)	2.75 ml	44 ml
Wash Buffer ** (Concentrate)	2 ml	25 ml
Elution Buffer	2.2 ml	30 ml
FABG Midi Column	2 pcs	25 pcs
Elution Tube (15 ml tube)	2 pcs	25 pcs
User Manual	1	1

Preparation of W1 Buffer and Wash Buffer for the first use:

Cat. No:	FABGK002-S (2 preps)	FABGK002 (25 preps)
* Ethanol volume for W1 Buffer	1 ml	16 ml
** Ethanol volume for Wash Buffer	8 ml	100 ml

Specification:

Principle: spin column (silica membrane)
Sample Size: up to 1.5 ml of fresh/frozen blood;
up to 6×10⁷ of cultured cells.

Column Capacity: 150 µg of DNA

Average DNA yield : $35\,\mu\text{g/ml}$ of whole blood


Handling Time: 1 hr Elution Volume: 1 ml

Required material to be provided by user

- •Pipettors and pipet tips
- •Centrifuge (should be capable up to 4,500 x g)
- Thermal incubator
- Oven (optional)
- •Ethanol (96~100%)
- Vortex

Important Notes:

- 1. Buffers provided in this system contain irritants. Wear gloves and lab coat when handling these buffers.
- 2. Preheat a thermal incubator to 60°C and 70°C before the operation.
- 3. Use a centrifuge with a swinging bucket rotor and a force of 4,500~6,000 x g for in all centrifugation steps.
- 4. Add ethanol (96~100%) to W1 Buffer and Wash Buffer before first use and store at room temperature.



General Protocol: For Whole Blood DNA Extraction

Please Read Important Notes Before Starting the Following Steps.

- 1. Transfer up to 1.5 ml sample (whole blood, buffy coat) to a 15 ml tube (not provided).

 -If lymphocytes sample is used, transfer $10^7 \sim 10^8$ cells to a 15 ml tube and adjust total volume to 1 ml with PBS.
- 2. Add 150 µl of Proteinase K to the sample and mix well by vortexing.
- 3. Add 1.5 ml of FABG Buffer to the sample and mix thoroughly by vortexing. Incubate the sample mixture at 60°C for 20 mins. During incubation, vortex briefly the tube 3 times and preheat Elution Buffer or ddH2O (0.5~1 ml per preparation) to 70°C.
 - -Do not add Proteinase K directly to FABG Buffer.
- 4. (Optional): If RNA-free genomic DNA is required, add 4 µl of 100 mg/ml RNase A (not provided) to the sample mixture and incubate at room temperature for 5 mins.
- 5. Add 1.5 ml of ethanol (96~100%) to the sample mixture and mix thoroughly by pulse-vortexing. If precipitate appears, break it by pipetting.
- 6. Place a FABG Midi Column to a 15 ml centrifuge tube (not provided). Transfer total sample mixture (ethanol added) carefully to the FABG Midi Column. Close the cap and centrifuge at 4,500~6,000 x g for 3 mins. Discard the flow-through and place the FABG Midi Column back to the 15 ml centrifuge tube.
- 7. Add 2 ml of W1 Buffer (ethanol added) to the FABG Midi Column. Close the cap and centrifuge at 4,500~ 6,000 x g for 3 mins. Discard the flow-through and place the FABG Midi Column back to the 15 ml centrifuge tube.
 - -Make sure that ethanol has been added into W1 Buffer at the first use.
- 8. Add 4.5 ml of Wash Buffer (ethanol added) to the FABG Midi Column. Close the cap and centrifuge at 4,500~6,000 x g for 10 mins. Discard the 15 ml centrifuge tube and the flow- through.
 - -Make sure that ethanol has been added into Wash Buffer at the first use.
 - -Avoid column tip touch the flow-through when transferring the FABG Midi Column.
 - -Note! 10 mins centrifugation is important for removing the residual of Wash Buffer from column membrane.
- 9. Transfer the FABG Midi Column to a new 15 ml centrifuge tube (Elution Tube, provided). Do not close the cap and stand the column at room temperature for 5 mins.
- 10. Add 0.5~1 ml of preheat Elution Buffer or ddH2O (pH 7.5-9.0) to the membrane center of FABG Midi Column. Stand the FABG Midi Column for 2 mins at room temperature.
 - -Important Step! Make sure that Elution Buffer is absorbed completely by column membrane.
- 11. Close the cap and centrifuge at 4,500~6,000 x g for 3 mins to elute DNA.

Protocol: For Cultured Cell DNA Extraction

Please Read Important Notes Before Starting the Following Steps.

- 1. Transfer up to 6×10^7 of cells to a 15 ml centrifuge tube (not provided). Centrifuge at 4,500 x g for 5 mins to pellet the cells.
 - If using adherent cells, trypsinize the cells before harvesting.
- 2. Resuspend the cells with 1.5 ml of PBS. Add 150 μ l of Proteinase K to the sample and mix well by vortexing.
- 3. Add 1.5 ml of FABG Buffer to the sample mixture and mix thoroughly by vortexing. Incubate the sample mixture at 60°C for 20 mins to lyse the sample. During incubation, invert the tube every 3~5 mins and preheat Elution Buffer or ddH2O (0.5~1 ml per preparation) to 70°C.
 - Do not add Proteinase K directly to FABG Buffer.
- 4. Follow the Whole Blood DNA Extraction protocol starting from step 4.